

## Dosage & Administration

### NEOVEMOX PLUS TAPE LONG ACTING HORSE WORMER & BOTICIDE GEL

*[Only relevant pack size data shown on label]*

The following table will act as a guide only:

Dose Table for 11.8g syringe:

Age	PONY		LIGHT HORSE		MEDIUM HORSE	
	Weight (kg)	Doses per syringe	Weight (kg)	Doses per syringe	Weight (kg)	Doses per syringe
3 months	60	9	125	4	155-170	3
6 months	95	6	185	3	230-265	2
Mature	200	2	400	1	500-575*	1

\*Note: One syringe is sufficient to treat a horse weighing 575 kg. For a heavy horse (weight exceeding 575 kg) it will be necessary to apply more than 1 syringe. For instance, a 725 kg horse would require 1¼ syringes to treat. To accurately determine the amount to administer with the second syringe, subtract 575 from the horse's weight and set the plunger accordingly. For a 725 kg horse, set the plunger on the second syringe to 150 kg (725 kg - 575 kg). Alternatively a 14.4 g syringe could be used to treat horses weighing up to 700 kg.

Dose Table for 14.4g syringe:

Age	PONY		LIGHT HORSE		MEDIUM HORSE	
	Weight (kg)	Doses per syringe	Weight (kg)	Doses per syringe	Weight (kg)	Doses per syringe
3 months	60	11	125	5	155-170	4
6 months	95	7	185	3	240-265	2
Mature	200	3	400	1	500-700*	1

\*Note: One syringe is sufficient to treat a horse weighing 700 kg. For a heavy horse (weight exceeding 700 kg) it will be necessary to apply more than 1 syringe. For instance, a 875 kg horse would require 1¼ syringes to treat. To accurately determine the amount to administer with the second syringe, subtract 700 from the horse's weight and set the plunger accordingly. For a 875 kg horse, set the plunger on the second syringe to 175 kg (875 kg - 700 kg).

## Dosage & Administration

### ADMINISTRATION

Prepare the syringe:

- Hold the syringe with the capped end pointing to the left so the weight measurements and the tick marks (small black lines);
- Turn the dial ring until the left side of the ring lines up with the weight of the animal;
- Remove the barrel cap and save it for later use if required.
- More than one syringe is needed to treat horses weighing more than 575 kg. [more than one syringe is needed to treat horses weighing more than 700 kg].

Administer the dose:

- Ensure the animal's mouth is empty;
- With your thumb on the plunger, insert the nozzle of the syringe into the horse's mouth at the space between front and rear teeth, and make sure the tip of the syringe is over the back of the tongue;
- Deliver the dose of gel onto the back of the tongue by gently pushing the plunger until it stops;
- Raise the horse's head once dosing is completed to ensure the gel is swallowed.

Surplus gel:

- Surplus gel can be stored or be used treat a second horse;
- If gel remains in the syringe once the full dose has been delivered, replace the barrel cap; Store the remaining gel below 30°C (room temperature) for later use, but do not use after the expiry date shown on the label.

If the first animal treated weighed less than 575 kg [700kg], there will be gel left in the syringe, which can be used to treat a second horse. To set the next dose, add the weight of the animal to be treated to the dose setting already on the syringe. For example if the syringe was first used to treat a foal weighing 100 kg, the green ring dial is set on 100 kg. To treat a horse of 400 kg next, move the green dial ring to the 500 kg marking.